

Emodren[®] cc

- A. Adaptor for silicone drainage tubes
- B. Multi-diameter connection
- C. White clamp
- D. Bellows 250 ml
- E. Red clamp (for discharge system)
- F. Collection bag



The Emodren Closed Circuit is a drainage system in continuous aspiration consisting of a 250 ml low-density PE bellows that generates a suction power of 100-110 mm Hg.

It has a no-return valve linked via a luer-lock to a 700 ml collection bag (it too fitted with an internal no-return valve) and a water-repellent air filter. The bag may be replaced one or more times with a spare one. The Emodren CC can be supplied in a base model without the drainage tube or in various versions with differing diameter drainage tubes (Ø 2.4 - 3.2 - 4.2 - 4.7 - 5.7 mm), with or without guide needles on them.

The Emodren CC makes it possible to drain off considerable quantities of fluid without disconnecting or opening the system, thus guaranteeing maximum safety for the operator, who may also choose whether to use a PVC or silicone drainage tube (the relevant connectors are in the kit provided). The Emodren CC is indicated for breast surgery, thyroidectomies, urological surgery, orthopaedic surgery, general surgery and whenever a washing/drainage procedure is required.

code	description	short description	qty
5003300	EMODREN CLOSED SYSTEM 250 ml	EMODREN C.C.	10
5143300	EMODREN C. C. with drainage tube Ø 2,4 mm Ch 7	EMODREN C.C. - TD14/22	10
5103300	EMODREN C.C. with drainage tube Ø 3,2 mm Ch 10	EMODREN C.C. - TD14/32	10
5113300	EMODREN C.C. with drainage tube Ø 4,2 mm Ch 13	EMODREN C.C. - TD12/42	10
5123300	EMODREN C.C. with drainage tube Ø 4,7 mm Ch 14	EMODREN C.C. - TD12/47	10
5133300	EMODREN C. C. with drainage tube Ø 5,7 mm Ch 17	EMODREN C.C. - TD12/57	10
5243300	EMODREN C.C. with drainage tube Ø 2,4 mm Ch 7 e ago guida	EMODREN C.C. - TD14/22 and needle	10
5203300	EMODREN C.C. with drainage tube Ø 3,2 mm Ch 10 e ago guida	EMODREN C.C. - TD14/32 and needle	10
5213300	EMODREN C.C. with drainage tube Ø 4,2 mm Ch 13 e ago guida	EMODREN C.C. - TD12/42 and needle	10
5223300	EMODREN C.C. with drainage tube Ø 4,7 mm Ch 14 e ago guida	EMODREN C.C. - TD12/47 and needle	10
5233300	EMODREN C.C. with drainage tube Ø 5,7 mm Ch 17 e ago guida	EMODREN C.C. - TD12/57 and needle	10
5053300	Compleat bag for EMODREN CLOSED SYSTEM	COLLECTION BAG	20