

Flexitip blade with detachable F.O.

- 4,0mm diameter F.O. bundle with more than 5.500 individual fiber strands
- Matt finished to reduce glare
- No chrome plating
- High quality fibers ensuring excellent light transmission, 7.000 lux with 2.5V Xenon/Halogen lamp, 12.500 lux with 3.5V Xenon/Halogen lamp
- Detachable fiber
- Movable tip facilitates difficult intubation
- Reduced likelihood of damage to teeth and soft tissues
- Autoclave up to 134°C e 4.000 standard cycles
- No adverse affect in using ETO, steam autoclaving or by soaking in Sporax or Cidex
- Can sustain plasma sterilization (Sterrad 100NX)
- The blades are compatible with handles ISO 7376



code	description	size
810-120-02	Flexitip F.O. blade	2
810-120-03	Flexitip F.O. blade	3
810-120-04	Flexitip F.O. blade	4

* On demand Set with handles LED or Xenon