

Macintosh blade single use

- Suitable of everyday hospital use
- Hygienic and secure
- Excellent light transmission, 4.000 lux with 2.5V Xenon/Halogen lamp, 10.000 lux with 3.5V Xenon/Halogen lamp
- Avoids the risk of infections
- The blades are compatible with handles ISO 7376



code	description	size
810-401-00	Macintosh blade disposable	0
810-401-01	Macintosh blade disposable	1
810-401-02	Macintosh blade disposable	2
810-401-03	Macintosh blade disposable	3
810-401-04	Macintosh blade disposable	4