

Macintosh F.O. blade with integrated fibers

- 4,5mm diameter F.O. with more than 6.500 individual fiber strands
- No chrome plating
- High quality fibers ensuring excellent light transmission, 8.000 lux with 2.5V Xenon/Halogen lamp, 14.000 lux with 3.5V Xenon/Halogen lamp
- Autoclave up to 134°C e 4.000 standard cycles
- No adverse affect in using ETO, steam autoclaving or by soaking in Sporax or Cidex
- Can sustain plasma sterilization (Sterrad 100NX)
- The blades are compatible with handles ISO 7376



code	description	size
810-001-00	Macintosh F.O. blade	0
810-001-01	Macintosh F.O. blade	1
810-001-02	Macintosh F.O. blade	2
810-001-03	Macintosh F.O. blade	3
810-001-04	Macintosh F.O. blade	4
810-001-05	Macintosh F.O. blade	5