

## Miller blade single use

- Suitable of everyday hospital use
- Hygienic and secure
- Excellent light transmission, 4.000 lux with 2.5V Xenon/Halogen lamp, 10.000 lux with 3.5V Xenon/Halogen lamp
- Avoids the risk of infections
- The blades are compatible with handles ISO 7376



code	description	size
810-411-00	Miller blade disposable	0
810-411-01	Miller blade disposable	1
810-411-02	Miller blade disposable	2
810-411-03	Miller blade disposable	3
810-411-04	Miller blade disposable	4
810-411-05	Miller blade disposable	00