

Miller blade with detachable F.O.

- 4,0mm diameter F.O. bundle with more than 5.500 individual fiber strands
- No chrome plating
- High quality fibers ensuring excellent light transmission, 7.000 lux with 2.5V Xenon/Halogen lamp, 12.500 lux with 3.5V Xenon/Halogen lamp
- Detachable fiber
- Autoclave up to 134°C e 4.000 standard cycles
- No adverse affect in using ETO, steam autoclaving or by soaking in Sporax or Cidex
- Can sustain plasma sterilization (Sterrad 100NX)
- The blades are compatible with handles ISO 7376



code	description	size
810-031-05	Miller F.O. blade	00
810-031-00	Miller F.O. blade	0
810-031-01	Miller F.O. blade	1
810-031-02	Miller F.O. blade	2
810-031-03	Miller F.O. blade	3
810-031-04	Miller F.O. blade	4