


Mini Emodren® 50/100

- A. Suction tube threaded connection
- B. Clamp
- C. Multi-diameter connection
- D. Drainage tube
- E. Bellows 50 ml
- F. Bellows 100 ml





With a relatively modest suction of 70-80 mm Hg, MINI EMODREN is most suited to surgery of small parts of the hand or foot, micro-surgery and plastic surgery.

The perforated drainage tube adapts perfectly to every anatomical profile, including fingers, thus preventing the formation of fibrous scar tissue in the operative wound and tendon sheaths.

It is supplied with a 7 cm suction tube, a manually detachable multi-diameter connection and a drainage tube TD 14/22 (Ø 2.4 mm).

The set can be provided with a drainage tube of a different diameter.

The two models (with 50 and 100 ml bellows) in the P version come packaged together with a disposable guide needle already assembled on the drainage tube.



code	description	short description	qty
10012592	MINI EMODREN 50 ml with drain tube Ø 2,4 mm Ch 7	EM 2/22 - TD 14/22	20
10012594	MINI EMODREN 50 ml with valve and with drain tube Ø 2,4 mm Ch 7	EM 2/22 CV - TD 14/22	10
10012593	MINI EMODREN P 50 ml with drain tube Ø 2,4 mm Ch 7 with guide-needle	EM 2/P - TD 14/22	10
10012612	MINI EMODREN P 50 ml with drain tube Ø 2,4 mm Ch 7 with short guide-needle	EM 2/P with short needle	10
CUSTOMMADE5	MINI EMODREN P 50 ml with drain tube Ø 3,2 mm Ch 10 with guide-needle	EM 2/P - TD 14/32	10
10012589	MINI EMODREN 100 ml with drain tube Ø 2,4 mm Ch 7	EM 2/100 - TD 14/22	10
CUSTOMMADE6	MINI EMODREN 100 ml with drain tube Ø 3,2 mm Ch 10	EM 2/100 - TD 14/32	10
CUSTOMMADE10	MINI EMODREN 100 ml with drain tube Ø 4,7 mm Ch 14	EM 2/100 - TD 12/47	10
10012590	MINI EMODREN P 100 ml with drain tube Ø 2,4 mm Ch 7 with guide-needle	EM 2/100 P - TD 14/22	10
CUSTOMMADE7	MINI EMODREN P 100 ml with drain tube Ø 3,2 mm Ch 10 with guide-needle	EM 2/100 P - TD14/32	10