



BLUESTAR PLUS Home Care with exhalation holes and security valve A.A.V

- Replaceable anatomic air cushion
- Wide range of fixing points
- Headgear mask in highly breathable material and with wide bearing surface for better forces distribution
- Operative openings for O2 enrichment and monitoring
- Port for the passage of the gastric tube, the introduction of the bronchoscope or for simple needs of the patient
- Safety valve with immediate opening in the event of air flow lack
- Exhalation holes for a correct evacuation of the exhaled breath
- Front support fully adjustable and removable

Bi-level CPAP VENTED

The presence of exhalation holes, safety valve Anti Asphyxia and replaceable air cushion, makes this mask the most versatile for ventilation of patients at home.

Each kit includes:

Mask with air cushion, transparent connector 22 mm M with security valve and exhalation holes, headgear mask, bulb pump

code	description	size
KM - 310 E	mask for N.I.V.	extra large
KM - 311 E	mask for N.I.V.	large
KM - 312 E	mask for N.I.V.	medium
KM - 313 E	mask for N.I.V.	small
KM - 314 E	mask for N.I.V.	pediatric