



BLUESTAR PLUS Intensive

- Replaceable anatomic air cushion
- Wide range of fixing points
- Headgear mask in highly breathable material and with wide bearing surface for better forces distribution
- Operative openings for O2 enrichment and monitoring
- Port for the passage of the gastric tube, the introduction of the bronchoscope or for simple needs of the patient
- Front support fully adjustable and removable

Bi-level CPAP NON-VENTED.

The design makes it widely versatile and the presence of 3 connectors in the Kit allows it to operate with different configurations of circuits. The connector with hole for exhalation allows to transform the mask from NON-VENTED to VENTED.

A port on the connector allows the O2 enrichment or monitoring.

Each kit includes:

Mask with air cushion, transparent connector 22 mm M, blue connector 22 mm F, exhalation valve, headgear mask, bulb pump.

code	description	size
KM - 310 B	mask for N.I.V.	extra large
KM - 311 B	mask for N.I.V.	large
KM - 312 B	mask for N.I.V.	medium
KM - 313 B	mask for N.I.V.	small
KM - 314 B	mask for N.I.V.	pediatric