



BLUESTAR PLUS Secure with Security valve A.A.V.

- Replaceable anatomic air cushion
- Wide range of fixing points
- Headgear mask in highly breathable material and with wide bearing surface for better forces distribution
- Operative openings for O₂ enrichment and monitoring
- Port for the passage of the gastric tube, the introduction of the bronchoscope or for simple needs of the patient
- Safety valve with immediate opening in the event of air flow lack
- Front support fully adjustable and removable

Bi-level CPAP Non-vented

The presence of a safety valve Anti-Asphyxia permits a reasonable airflow even in case of lack of air flow coming from the ventilator.

There is a swivel connector 360° blue color with a diameter of 22 mm F.

Each kit includes:

Mask with air cushion, blue connector 22 mm F with security valve, headgear mask, bulb pump.

code	description	size
KM - 310 V	mask for N.I.V.	extra large
KM - 311 V	mask for N.I.V.	large
KM - 312 V	mask for N.I.V.	medium
KM - 313 V	mask for N.I.V.	small
KM - 314 V	mask for N.I.V.	pediatric