

## Anaesthesia mask CRYSTAL™



- It fits the patient's face and optimizes the application of anesthetics
- PVC - DEHP and BFA free
- High transparency allows good visibility
- Safe and easy to use
- Single patient use
- Available in 8 sizes with color-coded fixing rings

- Each mask is packed individually with a description of the contents
- Appreciated in each country; guaranteed safety and comfort
- Vanilla scent

Up on request they can be delivered with a strawberry scent with the same code and adding the letter "S"



code	description	size	weight[gr]	dead space	color	connector	qty
KM-200	anaesthesia mask	extra large-7	59,5	310 ml	transparent	22F	50
KM-201	anaesthesia mask	large-6	55,5	250 ml	red	22F	50
KM-202	anaesthesia mask	medium-5	51,0	220 ml	yellow	22F	50
KM-203	anaesthesia mask	small-4	37,5	180 ml	green	22F	50
KM-204	anaesthesia mask	pediatric-3	27,0	85 ml	white	22F	50
KM-205	anaesthesia mask	infant-2	17,5	45 ml		15M	50
KM-206	anaesthesia mask	neonatal-1	14,5	25 ml		15M	50
KM-207	anaesthesia mask	pre-neonatal-0	11,0	20 ml		15M	50