

Manual resuscitator

With overpressure control pop off valve and mask

- PVC and PC, DEHP and BFA free resuscitator bag
- Latex Free
- Adults 1700ml volume
- Pediatric 700ml volume
- Infant 300ml volume
- Single patient use



Kit with: resuscitator bag, oxygen tube, reservoir, Crystal mask and overpressure "pop off" control valve

code	description	size	with mask	qty
KM-806A	manual resuscitator	adult	KM-202 medium	10
KM-807A	manual resuscitator	pediatric	KM-204 pediatric	10
KM-808A	manual resuscitator	infant	KM-205 infant	10