

Manual resuscitator

With overpressure control pop off valve and PEEP valve adapter

- PVC and PC, DEHP and BFA free resuscitator bag
- Latex Free
- Adults 1700ml volume
- Pediatric 700ml volume
- Infant 300ml volume
- Single patient use



Kit with: resuscitator bag, oxygen tube, reservoir, without mask, overpressure "pop off" control valve and PEEP valve adapter

code	description	size	qty
KM-8061	manual resuscitator	adult	10
KM-8071	manual resuscitator	pediatric	10
KM-8081	manual resuscitator	infant	10