

## Mask for CPAP Zephir



- Safety valve
- Inflating cushion with inflation valve. It fits perfectly ensuring a high hold and avoiding air leak
- Provided with port for oxygen and monitoring
- In PVC -DEHP and BPA free high transparency allows good visibility
- Latex Free
- Safe and easy to use
- Provided with two ports with valve
- Single patient use
- Available in 3 sizes

code	description	size
KM-301A	mask for CPAP	large
KM-301DA	mask for CPAP	large deep
KM-303A	mask for CPAP	small