

Skin Stapler

- Simple and intuitive use
- Highly biocompatible stainless steel staples
- Secure due to the return valve
- Staple automatic release system
- Non-stop recharge
- Less bleeding and less risk of cross-contamination



The disposable Skin Stapler by MeHoS makes suturing a quick, easy and safe process. It helps the surgeon considerably and gives to the patient a quick and pain-less healing.

Thanks to the automatic release system of the staple the surgeon does not need to use any pressure on the patient's skin, pressure that would make the removal of the staples painful.

The staples are made in medical grade stainless steel and they can stay in place for a long time, thanks to their biocompatibility, giving the patient a quicker recovery.

code	description	size	∅ mm	clip size after the application	qty
F-10WJ	Skin stapler with 10 staplers	wide	0.6	width 7.2 mm, height 4.9 mm	5
F-35WE	Skin stapler with 35 staplers	wide	0.6	width 7.2 mm, height 4.9 mm	5
F-35RE	Skin stapler with 35 staplers	regular	0.5	width 5.9 mm, height 3.9 mm	5
F-55WE	Skin stapler with 55 staplers	wide	0.6	width 7.2 mm, height 4.9 mm	5
F-55RE	Skin stapler with 55 staplers	regular	0.5	width 5.9 mm, height 3.9 mm	5



Remover

- Latex Free
- Disposable
- Sterile

code	description	qty
C-WI	Remover	20